



AB Diagnostic Systems

REF

DK.064.01.8



96

IVD

For *In vitro* Diagnostic Use

INSTRUCTIONS FOR USE
abia Chlamydia Ab IgA
Enzyme immunoassay for the detection
of IgA antibodies to *Chlamydia trachomatis*
in human serum or plasma

This Package Insert provides information for Professional Use of the kit.

The kit contains sufficient reagents for 96 assays (one breakable plate) including controls; the kit is intended for manual testing with a possibility of fractional (one strip) use of the kit or use of the kit on open type automated analyzer for enzyme immunoassay.

I. INTENDED USE

The abia Chlamydia Ab IgA is an enzyme immunoassay for the detection of IgA antibodies specific to *C. trachomatis* in human serum (plasma). The kit is used as an aid in the diagnosis of *C. trachomatis* specific infection.

The abia Chlamydia Ab IgA is intended to be run and interpreted in conjunction with the abia Chlamydia Ab IgG kit and the abia Chlamydia Ab IgM.

This kit is for diagnostic use by a trained laboratory professional and will not be sold to the general public. All the reagents are for professional *in vitro* diagnostic use only.

The results of this or any other diagnostic assay should be used and interpreted only in the context of the overall clinical picture.

II. INTRODUCTION

Chlamydia is a gram negative obligate intracellular bacteria that causes acute and chronic disease in mammalian and avian species. The genus Chlamydia is comprised of four species: *C. trachomatis*, *C. pneumoniae*, *C. psittaci* and *C. pecorum*. *C. trachomatis* is divided into 15 serovars. Serovars A, B, Ba and C are agents of trachoma, the leading cause of preventable blindness, endemic in third world countries. Serovars L1-L3 are the agents of lymphogranuloma venereum. Serovars D-K are the common cause of sexually transmitted genital infection worldwide: cervicitis, endometritis/salpingitis in females and urethritis in both males and females.

Serological cross reactions occur between the three different species of Chlamydia: *C. trachomatis*, *C. pneumoniae* and *C. psittaci*. Most of the serological diagnostic assays for Chlamydia use either purified elementary bodies: microimmunofluorescence (MIF) and EIA tests, lipopolysaccharide (LPS), or purified major outer membrane protein (MOMP) as antigens. Genus specific epitopes are present in all the above antigens, therefore, low species specificity is observed. Moreover, a large proportion of the population has been exposed to *C. pneumoniae* (with no clinical signs), the prevalence of anti-Chlamydia antibodies is very high. Therefore, the differentiation between *C. pneumoniae* and *C. trachomatis* specific antibodies using conventional serological screening tests (MIF, EIA, EIA etc.) is insufficient.

C. trachomatis species specific epitopes, derived from MOMP, are used in an Enzyme Linked Immunosorbent Assay. The test excludes cross-species reactive epitopes and enables more accurate and more specific determination of *C. trachomatis* IgG, IgA and IgM antibodies.

III. PRINCIPLE OF THE TEST

The abia Chlamydia Ab IgA is an indirect two-step immunoassay for the detection of IgA antibodies to *C. trachomatis*. The plates are coated with *C. trachomatis* specific proteins. Serum to be tested is diluted and incubated with the precoated plate. In this step *C. trachomatis* specific antibodies are bound to the immobilized *C. trachomatis* specific proteins. Non specific antibodies are removed by washing. Anti-human IgA conjugated with horseradish peroxidase (HRP) is added and incubated. In this step the HRP-conjugate is bound to the prebound antigen-antibody complex. Unbound conjugate is removed by washing. The presence of bound enzyme indicating the presence in the specimen of specific antibodies is revealed by a color change in substrate/chromogen (TMB) mixture.

IV. CONTENT OF THE KIT abia Chlamydia Ab IgA

Table 1

LABEL	NATURE OF THE REAGENTS	PRESENTATION
Chlamydia TR-Ag Coated Strips	Polystyrene stripped 96-well plate (breakable wells) coated with a mix of recombinant proteins, which represent the recombinant analogs of the MOMP. Store at 2-8 °C until expiration date.	1 plate
Conjugate (concentrated 21-fold)	Antibodies against human IgA, conjugated with HRP enzyme with addition of bovine serum albumin (1.73%). Glycerol based solution. Transparent or slightly opalescent colorless or light yellow liquid. Preserving agent: 0.04% ProClin 300, 0.0009% gentamycin sulfate. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 0.75 ml
Positive Control, Inactivated	Inactivated human serum, containing IgA antibodies to <i>Chlamydia trachomatis</i> . The serum does not contain HBsAg, antigen p24 HIV-1, HIV-1,2 and HCV antibodies. Transparent or slightly opalescent crimson-red liquid. Preserving agents: 0.10% ProClin. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 1.5 ml
Negative Control, Inactivated	Inactivated human serum, not containing IgA antibodies to <i>Chlamydia trachomatis</i> . The serum does not contain HBsAg, antigen p24 HIV-1, HIV-1,2 and HCV antibodies. Transparent or slightly opalescent green liquid. Preserving agents: 0.04% ProClin 300, 0.19% sodium azide. 0.001% gentamycin sulfate. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 3.0 ml
Sample Diluent	Sample buffer that is used to dilute samples before analysis. Transparent violet-blue liquid with possible sedimentation. Preserving agent: 0.01% thimerosal. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 12.5 ml
Conjugate Diluent	Buffer that is used to dilute Conjugate concentrated before analysis. Transparent or slightly opalescent yellow liquid with possible sedimentation. Preserving agent: 0.01% thimerosal. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 13.5 ml
Washing Solution (concentrated 25-fold)	Phosphate-saline solution (pH 7.4-7.7). Transparent or slightly opalescent liquid, colorless or pale yellow, sediment may form that dissolves at 35-39 °C and shaking. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 50.0 ml

Substrate Buffer	Citric acid (0.64%) solution, pH 4.1-4.3, containing H ₂ O ₂ (0.008%). Transparent colorless liquid. Preserving agent: 0.04% ProClin 300. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 15.0 ml
TMB (concentrated 11-fold)	Solution containing 3,3',5,5'-Tetramethylbenzidine (TMB) (0.4%) and dimethyl sulfoxide (DMSO) (85.60%). Transparent colorless liquid. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 1.5 ml
Stopping Reagent	Sulfuric acid solution (H ₂ SO ₄) 0.75M. Transparent colorless liquid. Store at 2-8 °C until expiration date in a tightly sealed vial.	1 vial 25.0 ml
Protective films for EIA plates		1
Disposable tips		16
Disposable plastic dishes for liquid reagents		2
Polyethylene bag with a Zip-Lock		1

V. PRECAUTIONS

The reliability of the results depends on correct implementation of the following requirements:

- The temperature in the lab should be 18-24 °C.
- Inspect the contents of the box: check the vials and labels integrity. In case of label loss or labels/vials damage, vials should be disposed and **kit cannot be used**.
- Do not use expired reagents.
- Do not mix reagents from different lots within a given test run.
- Carefully reconstitute the reagents avoiding any contamination.
- Do not carry out the test in the presence of reactive vapors (acid, alkaline, aldehyde vapors) or dust that could alter the enzyme activity of the conjugates.
- Use glassware thoroughly washed and rinsed with deionized water or preferably, disposable material.
- Do not allow the microplate to dry between the end of the washing operation and the reagent distribution.
- The enzyme reaction is very sensitive to metal ions. Consequently, do not allow any metal element to come into contact with the various conjugate or substrate solutions.
- Use a new distribution tip for each sample.
- Do not reuse protective films for EIA plates.
- Well washing is a critical step in this procedure: respect the recommended number of washing cycles and make sure that all wells are completely filled and then completely emptied. Incorrect washing may lead to inaccurate results.
- Never use the same container to distribute conjugate and other solutions.
- Check the pipettes and other equipment for accuracy and correct operation.
- Do not change the assay's procedure.
- Use distilled or deionized water.
- Avoid exposure of the reagents to excessive heat or sunlight during storage and incubation.
- Once the assay has been started, all subsequent steps should be performed without interruption.

VI. HEALTH AND SAFETY INSTRUCTIONS

- All reagents included in the kit are intended for “in vitro diagnostic use”.
- Human origin material used in the preparation of Negative Control and Positive Control has been tested and found negative for HBsAg, antigen p24 HIV-1, antibodies to hepatitis C virus and antibodies to human immunodeficiency virus (HIV-1 and HIV-2).
- Certain reagents contain biological material of animal origin.
- Because no known test method can offer complete assurance that infectious agents are absent, handle reagents and patients samples as if capable of transmitting infectious disease.
- Do not eat, drink, smoke or apply cosmetics where immunodiagnostic materials are being handled.
- Any equipment directly in contact with specimens and reagents as well as washing solutions should be considered as contaminated products and treated as such.
- Wear lab coats and disposable gloves when handling reagents and samples and thoroughly wash your hands after handling them.
- Avoid spilling samples or solutions containing samples.
- Avoid any contact of Substrate Buffer, TMB and Stopping Reagent with the skin and mucosa.
- Provide adequate ventilation.
- All materials contacted with specimens or reagents, including liquid and solid wastes, should be inactivated by validated procedures (autoclaving or chemical treatment) and disposed in accordance with applicable local law regulations.



Warning!

Positive Control contains ProClin 300.
H317: May cause an allergic skin reaction.
P261: Avoid breathing vapors.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352: IF ON SKIN: Wash with plenty of water.
P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.



Warning!

Negative Control contains sodium azide.
H312: Harmful in contact with skin.
P280: Wear protective gloves/protective clothing/ eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of water. Immediately call a poison center/doctor.
P312: Call a POISON CENTER or doctor/physician if you feel unwell.



Danger!

Stopping Reagent contains 0.75M sulfuric acid.
H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P303 + P361 + P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor/ physician.

VII. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED WITH THE KIT:

- Distilled or deionized water.
- Automatic or semiautomatic, adjustable or preset single-channel and multi-channel pipettes with a changeable volume for a set of liquids.
- Disposable pipette tips.
- Microplate incubator at $(37.0 \pm 1.0) ^\circ\text{C}$.
- Automatic microplate washer.
- Microplate reader equipped with 450 nm or with 450 and 620-680 nm filters.
- Open type automated analyzer with 450 nm filter or with 450 and 620-680 nm filters (for automated procedure).
- Laboratory clock.

VIII. COLLECTION AND HANDLING OF SPECIMENS

Blood samples should be collected according to the current practices. Serum, plasma may be used. Separate serum or plasma from blood cells as soon as possible to avoid any haemolysis. Extensive haemolysis may affect test performance. Specimens with observable particulate matter should be clarified by centrifugation prior to testing. Suspended fibrin particles or aggregates may yield falsely positive results. Do not heat the samples.

Store/transport the samples in accordance with the current regulatory documentation. If samples are to be stored/transported for a longer period of time, they must be frozen at or below $-20 ^\circ\text{C}$. Avoid repeated freeze/thaw cycles. Samples that have been frozen and defrosted more than 1 time cannot be used. Samples with expressed haemolysis, hyperlipidemia must not be analyzed.

IX. PREPARATION OF THE REAGENTS

1. Ready to use reagents:

- **Chlamydia TR-Ag Coated Strips.** Each plate containing 12 strips is wrapped in a sealed foil-lined bag. Open the bag and remove the plate. Select the number of coated strips required for the assay. Unused strips should be placed back into the bag. After the bag has been opened, the strips are stable within 1 month after opening when stored at $2-8 ^\circ\text{C}$, provided that the foil pack is resealed in polyethylene bag with a Zip-Lock. The silica gel bag should not be removed from the foil packaging.
- **Positive Control;**
- **Negative Control;**
- **Sample Diluent;**
- **Stopping Reagent (0.75M).**

2. Reagents to prepare:

- **Working Washing Solution.** Thoroughly mix the contents of the bottle with concentrated Washing Solution (concentrated 25-fold). Dilute the required volume of concentrated Washing Solution with the corresponding volume of distilled or deionized water prior to use (See Table 2). Mix the solution thoroughly.
- **Working Solution of Conjugate** must be prepared before usage. Thoroughly mix Conjugate concentrate. To make Working Solution of Conjugate, take required amount of concentrate and mix with Conjugate Diluent (See Table 2) in a separate vial. Mix thoroughly the content of the vial avoiding foaming (Do not use intensive mixing!).
- **Substrate Mixture** must be prepared before usage. Take required amount of TMB (concentrated 11-fold) to a separate vial and add required amount of Substrate Buffer (See Table 2). Mix thoroughly.

Substrate mixture should be colorless!

The volume of reagents required for the certain number of strips or plate is provided in the table below:

Table 2

Reagent preparation

Number of strips to be used	Working Washing Solution		Working Solution of Conjugate		Substrate Mixture	
	Washing Solution (concentrated 25-fold) (ml)	Distilled or deionized water (ml)	Conjugate (concentrated 21-fold) (ml)	Conjugate Diluent (ml)	TMB (concentrated 11-fold) (ml)	Substrate Buffer (ml)
1	4.0	96.0	0.05	1.0	0.1	1.0
2	8.0	192.0	0.10	2.0	0.2	2.0
3	12.0	288.0	0.15	3.0	0.3	3.0
4	16.0	384.0	0.20	4.0	0.4	4.0
5	20.0	480.0	0.25	5.0	0.5	5.0
6	24.0	576.0	0.30	6.0	0.6	6.0
7	28.0	672.0	0.35	7.0	0.7	7.0
8	32.0	768.0	0.40	8.0	0.8	8.0
9	36.0	864.0	0.45	9.0	0.9	9.0
10	40.0	960.0	0.50	10.0	1.0	10.0
11	44.0	1056.0	0.55	11.0	1.1	11.0
12	50.0	1200.0	0.65	13.0	1.2	12.0

X. TEST PROCEDURE

Note: Before use, allow reagents to reach room temperature (18-24 °C) for 30 min.

Step	The assay procedure
1	Wash the plate with Working Washing Solution twice before assay run. Carefully fill all wells with at least 380 µl of Working Washing Solution with multichannel pipette or automatic microplate washer, wait for 40 seconds and remove Working Washing Solution from the wells. Do not leave any fluid in the wells.
2	Add 100 µl of Positive Control, Negative Control into the wells of Coated Strips. <u>1 strip</u> – Positive Control to 1 well, Negative Control to 1 well; <u>2 strips</u> – Positive Control to 1 well, Negative Control to 2 wells; <u>3 strips and more</u> – Positive Control to 2 wells, Negative Control to 2 wells.

3	Add 90 µl of Sample Diluent and 10 µl of the sera samples to the rest of the wells (serum dilution ratio is 1:10). Carefully mix fluid in wells by gentle pipetting. Violet-blue color should change to blue-green. If you do not observe change of the color then test results may be false, or there is no serum added to the well.
4	Cover the strips with a protective film. Incubate 30 min in a microplate incubator at (37.0 ± 1.0) °C.
5	Remove liquid from the wells, wash the plate 4 times as described in step 1.
6	Add 100 µl of Working Solution of Conjugate into all wells.
7	Cover the plate with a protective film. Incubate for 30 min in a microplate incubator at (37.0 ± 1.0) °C.
8	Remove liquid from the wells, wash the plate 5 times as described in step 1.
9	Add 100 µl of Substrate Mixture into all wells.
10	Incubate at 18-24 °C for a 10 min in a dark place.
11	Add 50 µl of Stopping Reagent into each well to stop the reaction and read the optical density at 450/620-680 nm using a microplate reader. Reading the absorbance at 450 nm only is possible.

Scheme of the assay is represented in Annex.

Automated analyzers

For automated test procedure, it is advisable to use protocol submitted by the manufacturer. When creating the protocol independently, follow the procedure specified in section X TEST PROCEDURE, and comply with the requirements provided in sections V. PRECAUTIONS.

When preparing working solutions of reagents for the automated test procedure, dead volume of vials or containers used to place the solutions onboard should be taken into account.

Validated test protocols and dilution tables of working solutions for different models of EIA analyzers can be obtained upon request from the manufacturer (see section XV).

XI. RESULTS

The presence or absence of antibodies against *Chlamydia trachomatis* is determined by the ratio of the OD of each sample to the calculated Cut-Off value.

Test Validation

For the test to be valid the following criteria must be met. If these criteria are not met the test should be considered invalid and should be repeated.

1. **Positive Control:** The absorbance value should not be less than 0.800.
2. **Negative Control:** The absorbance value should not be more than 0.200.

Calculate Cut-Off value as:

$$\text{Cut-Off} = \text{average OD value of Negative Control} + A \quad (A = 0.220),$$

where **A**– is a coefficient defined by manufacturer during statistical processing for each lot.

Interpretation of Result

Sample is positive, if the OD value is \geq Cut-Off.

Sample is negative, if the OD value is $<$ Cut-Off.

Determination of Antibody Titer

The abia Chlamydia Ab IgA can be used to reveal titer of anti-CHLAMYDIA TR IgA antibodies in human serum (plasma). Positivity Index (PI) is used to reveal titer:

$$PI = OD (\text{sample}) / \text{Cut-Off}$$

Then Table “Correlation of PI and anti-CHLAMYDIA TR IgA antibodies titer” is used to determine antibodies titer.

Correlation of PI and anti-CHLAMYDIA TR IgA antibodies titer

Positivity Index	Serum titer
1-2.5	1/10
2.6-5.0	1/20
5.1-8.0	1/40
8.1-11.9	1/80
≥ 12.0	$\geq 1/160$

XII. PERFORMANCE CHARACTERISTICS OF abia Chlamydia Ab IgA

1. Concordance

The abia Chlamydia Ab IgA was compared to another commercially available immunoassay as a reference test. The total number of specimens was 215. Concordance = 91.4%.

XIII. LIMITS OF THE TEST

1. No single serological test should be used for a final diagnosis. All clinical and laboratory data should be taken into account.

2. Samples obtained too early during primary infection may not contain detectable antibodies. If Chlamydial infection is suspected, a second sample should be obtained 14-21 days later and tested in parallel with the original sample.

XIV. CONDITIONS OF STORAGE AND TRANSPORTATION

- **Expiry date is indicated on the packaging.** Storage and transportation conditions for the kit, conditions and terms of storage for working solutions and unused reagents are specified in Table 3.
- Transportation should be done by covered transport at specified temperature in accordance with established transportation regulations. Kits transported at improper temperature cannot be used.
- Kits stored improperly cannot be used.

Table 3

1	Storage conditions		
	Keep in a dark dry place at 2-8 °C. Freezing is prohibited.		
2	Transportation conditions		
	at 2-8 °C		
	at 9-20 °C	not more than during ten (10) days	
3	Conditions and terms of storage for working solutions		
	Keep in a dark dry place and in a chemically neutral vial		
	Working Washing Solution	at 2-8 °C	For up to 3 days
	Working Solution of Conjugate	at 18-24 °C in a dark place	For up to 12 hours
	Substrate Mixture	at 18-24 °C in a dark place	For up to 10 hours
4	Conditions and terms of storage of unused reagents after opening		
	Keep in a dark dry place at 2-8 °C.		
	Chlamydia TR-Ag Coated Strips	Place the unused strips/wells back into the bag, reseal the foil-lined package in Zip-Lock plastic bag. Do not remove desiccant.	For 1 month
	Positive Control, Negative Control, Conjugate, Sample Diluent, Conjugate Diluent, Washing Solution, Substrate Buffer, TMB, Stopping Reagent	Close the vials tightly with screw caps and store them in the manufacturer's package.	Until the kit expiration date

XV. GUARANTEE

- Manufacturer guarantees conformity of the product to the requirements of regulatory and technical documentation.
- Quality and safety of the kit is guaranteed within established shelf life.
- Please contact Manufacturer if you have any questions.



AB Diagnostic Systems GmbH

Sportfliegerstraße 4, Berlin, 12487, Germany

Tel. +49 30 208987160, Fax: +49 30 208987199

E-Mail: info@ab-ds.de, www.ab-ds.de

XVI. REFERENCES

1. Immunoreactivity studies on synthetic peptides deriving from variable domain IV of *Chlamydia trachomatis* major outer membrane protein. S. Klimashevskaya, T. Ulanova, A. Burkov, A. Obriadina. 17th European Congress of Clinical Microbiology and Infectious Disease. - Munich, Germany, 2007. P. 1140.

2. Antigen specific serum antibody response to *Chlamydia trachomatis* in patients with acute pelvic inflammatory disease. A. Miettinen, P.K. Heinonen, K. Teisala, R. Punnonen, J. Paavonen. J. Clin. Pathol. 1900; 43:758-761.

3. Antibody recognition of a neutralization epitope on the major outer membrane protein of *Chlamydia trachomatis*. Zhong G, Berry J, Brunham RC. Infect. Immun. 1994 May; 62(5):1576-83.
















4. Characterization of the disulfide bonds and free cysteine residues of the *Chlamydia trachomatis* mouse pneumonitis major outer membrane protein. Yen T.Y., Pal S., de la Maza L.M. Biochemistry. 2005 Apr 26; 44(16):6250-6.

5. Functional and structural mapping of *Chlamydia trachomatis* species-specific major outer membrane protein epitopes by use of neutralizing monoclonal antibodies. Peterson E.M., Cheng X., Markoff B.A., Fielder T.J., de la Maza L.M. Infect. Immun. 1991 Nov; 59(11):4147-53.

6. Comparison of three commercially available peptide-based immunoglobulin G (IgG) and IgA assays to microimmunofluorescence assay for detection of *Chlamydia trachomatis* antibodies. Morre S.A., Munk C., Persson K., Kruger-Kjaer S., van Dijk R., Meijer C.J., van Den Brule A.J. J. Clin. Microbiol. 2002 Feb; 40(2):584-7.

7. *Chlamydia trachomatis* serology: diagnostic value of outer membrane protein 2 compared with that of other antigens. Bas S, Muzzin P, Vischer T.L. J. Clin. Microbiol. 2001 Nov; 39(11):4082-5.

XVII. EXPLANATION OF SYMBOLS

	Manufacturer		Storage temperature limitation
	Date of manufacture CCYY-MM		Consult Instruction for use
	Expiry date CCYY-MM-DD		For in vitro diagnostic use
	Batch code		Sufficient for
	Catalog number		Symbol "exclamation mark"
	Fragile, handle with care	Warning!	Signal word
	Keep away from sunlight		Symbol "corrosion"
	Keep dry	Danger!	Signal word
	Top		

Scheme of the assay

1	Wash the plate	Working Washing Solution, not less than 380 µl, 2 times
2	Add	100 µl of Positive Control, Negative Control (wells of the Coated Strips) <u>1 strip</u> – Positive Control to 1 well, Negative Control to 1 well; <u>2 strips</u> – Positive Control to 1 well, Negative Control to 2 wells; <u>3 strips and more</u> – Positive Control to 2 wells, Negative Control to 2 wells.
3	Add	90 µl of Sample Diluent and 10 µl of samples
4	Incubate	30 min, (37.0 ± 1.0) °C, microplate incubator
5	Wash the plate	Working Washing Solution, not less than 380 µl, 4 times
6	Add	100 µl of Working Solution of Conjugate
7	Incubate	30 min, (37.0 ± 1.0) °C, microplate incubator
8	Wash the plate	Working Washing Solution, not less than 380 µl, 5 times
9	Add	100 µl of Substrate Mixture
10	Incubate	10 min, 18-24 °C in the dark place
11	Add	50 µl of Stopping Reagent
12	Read the optical density	450 nm/620-680 nm or 450 nm