

# abia Progesterone




REF DK.039.01.3

IVD



Note: Changes highlighted ★



# abia

## Intended use

Abia Progesterone enzyme immunoassay for the quantitative determination of progesterone concentration in human serum.

For professional use only.

## Clinical value

Progesterone is a female sex hormone which, in conjunction with estrogens, regulates the accessory organs during the menstrual cycle and it is particularly important in preparing the endometrium for the implantation of the blastocyte and in maintaining pregnancy.

In non-pregnant women progesterone is mainly secreted by the corpus luteum whereas in pregnancy the placenta becomes the major source. Minor sources are the adrenal cortex for both sexes and the testes for males. Progesterone circulates in blood mainly bound to corticosteroid binding globulin (CBG), sex hormone binding globulin (SHBG) and albumin. Only 2 - 10 % of the total concentration circulates as free hormone. Blood progesterone concentrations vary widely according to the phases of menstrual cycle. The maximal levels are achieved 4 - 7 days after ovulation and remain elevated for 4 - 6 additional days prior to falling to the preovulatory levels 24 hours before the onset of menstruation.

Since the rise and fall of progesterone parallel the activity of ovarian follicle and corpus luteum, measurements of serum progesterone are clinically used to confirm ovulation and normal function of the corpus luteum in non-pregnant women. If ovulation does not occur the corpus luteum is not formed and no cyclical rise of progesterone in blood is observed.

Abnormal progesterone secretion has been implicated in premenstrual tension, irregular shedding of endometrium, dysmenorrhoea, and luteal insufficiency. Progesterone concentration can vary not only from subject to subject but also in the same person from day to day or even

from hour to hour. Consequently, in gynecological disorders or abnormal pregnancies serial measurements rather than single ones are recommended for a proper interpretation of results. During pregnancy progesterone is widely produced by placenta.

## Principle of the test

Abia Progesterone is a one-step immunoassay, based on the principle of the competitive method.

Progesterone present in the sample and the labeled enzyme-progesterone in the conjugate compete for binding to the capture antibody on the anti-progesterone coated microplate. The enzyme activity in the antibody-bound fraction is inversely proportional to the native progesterone concentration.

The unbound components are removed by washing. After addition of the solution containing TMB and hydrogen peroxide, the wells with bound conjugate develop a blue color which is converted to yellow after the reaction has been stopped with sulphuric acid.

The color intensity is inversely proportional to the concentration of the progesterone in the specimen and can be read at 450 nm.

## Kit contents

	S	
Progesterone Ab coated plate	1	polystyrene plate 12 × breakable 8-well strips coated with monoclonal anti-progesterone antibodies
Conjugate	1 × 12 ml	ready to use; HRP-labeled progesterone; transparent or slightly opalescent pink liquid
Calibrator 0	1 × 0.5 ml	human serum not containing progesterone; pale yellow liquid
Calibrator 1	1 × 0.5 ml	human serum containing progesterone in concentration approx. 1.5 nmol/l; pale yellow liquid
Calibrator 2	1 × 0.5 ml	human serum containing progesterone in concentration approx. 10 nmol/l; pale yellow liquid
Calibrator 3	1 × 0.5 ml	human serum containing progesterone in concentration approx. 25 nmol/l; pale yellow liquid
Calibrator 4	1 × 0.5 ml	human serum containing progesterone in concentration approx. 50 nmol/l; pale yellow liquid
Calibrator 5	1 × 0.5 ml	human serum containing progesterone in concentration approx. 100 nmol/l; pale yellow liquid
Control serum	1 × 0.5 ml	human serum based control containing progesterone; pale yellow liquid
Washing solution (concentrated 25-fold)	1 × 50 ml	phosphate saline buffer; colorless or pale yellow liquid
TMB/substrate solution	1 × 12 ml	ready to use; citric acid buffer containing TMB and H <sub>2</sub> O <sub>2</sub> ; colorless liquid
Stopping reagent 0.2M H <sub>2</sub> SO <sub>4</sub>	1 × 25 ml	ready to use; 0.20 mol/l sulphuric acid solution; colorless liquid
Protective film	1	
Plastic dish	2	
Plastic zip-lock bag	1	

Exact concentration levels for calibrators and control serum are given on the labels and certificates of analysis on a lot specific basis. For conventional units: nmol/l x 0.314 = ng/ml.

All components are stable until expiration date of the kit when stored at 2–8 °C in a tightly sealed package. Expiration date is indicated on the package.

Once opened, the components should be used within 2 months. Concentration of preserving agents: <=0.1 %.

## Materials and equipment required but not provided

- purified water
- automatic or semiautomatic, adjustable or preset pipettes or multipipettes
- disposable pipette tips
- automatic microplate washer
- microplate reader equipped with 450 nm filter

## Safety notes

- human origin material used in the preparation of the calibrators and control serum has been tested and found non reactive for hepatitis B surface antigen (HBsAg), antigen p24 HIV-1, antibodies to hepatitis C virus and antibodies to human immunodeficiency virus (HIV-1 and HIV-2)
- as no known test method can offer complete assurance that infectious agents are absent, reagents and samples should be handled as if capable of transmitting infectious disease; any equipment directly in contact with samples and reagents should be considered as contaminated
- do not eat, drink, smoke or apply cosmetics in the laboratory
- do not pipette by mouth
- avoid any contact of the reagents and samples with the skin and mucosa; wear lab coats and disposable gloves when handling them; thoroughly wash your hands after work
- avoid spilling samples or solutions containing samples. Wipe spills immediately and decontaminate affected surfaces
- all materials contacted with specimens or reagents, including liquid and solid waste, should be inactivated by validated procedures (autoclaving or chemical treatment) and disposed in accordance with applicable local law regulations

## Precautions

- do not use reagents without label or with damaged label/package
- do not use expired reagents
- do not change the assay procedure; perform all subsequent steps without interruption
- do not mix reagents from different lots
- do not mix the caps of vials
- do not run the EIA test in the presence of reactive vapours (acid, alkaline, aldehyde), dust or metals
- do not let the wells dry once the assay has been started
- do not use the same container and tips for different liquid components of the kit and samples
- do not reuse the coated plates
- do not reuse the removed protective film
- do not expose the reagents to excessive heat or sunlight during storage and test procedure
- do not freeze the reagents

## Collection and handling of specimens

- collect blood specimens according to the current practices
- use serum for testing; performances of the test have not been evaluated on other biological fluids
- separate the clot or red cells from serum as soon as possible to avoid any haemolysis
- do not use sera preserved with sodium azide, thiomersal
- do not use contaminated, hyperlipaemic and hyperhaemolysed specimens
- the samples with hyperproteinaemia and hyperbilirubinaemia were not specially tested
- before testing samples with observable particulate matter should be clarified by centrifugation
- suspended fibrin particles or aggregates may yield reactive results
- do not heat the samples
- samples can be stored at 2–8 °C within 72 hours or deep-frozen at -20 °C
- no more than one freeze/thaw cycle is allowed

## Procedural notes

- before use wait 30 minutes for the reagents to stabilize to room temperature (20–25 °C)
- check appearance of the reagents
- lost vacuum in the bag of the coated plate will not affect the performance of the test
- check the pipettes and other equipment for accuracy and correct operation
- the washing procedure is a critical step; for the detailed washer settings see section “Washing procedure”
- for the description of test procedure with the automated analyzers see section “Automated analyzers”

## Washing procedure

Please contact your representative for protocols for recommended washers and procedures. In general the following protocol is recommended:

- flow-through washing with a volume not less than 300 µl per well is used
- repeat 5 times
- do not allow the wells to become dry during the assay procedure
- ensure that no liquid is left in the well (use double aspiration in the final step where possible)
- avoid to tap out the plate
- residual volume lower than 10 µl is not critical for following steps of the test procedure
- when using a microplate washer clean the wash head frequently to prevent contamination

## Preparation of reagents

Number of strips to be used	1	2	3	4	5	6	7	8	9	10	11	12
<b>Working washing solution:</b> mix the reagents thoroughly by inversion <b>Stability:</b> 14 days at 18–24 °C or 28 days at 2–8 °C												
Washing solution (concentrated 25-fold), ml	3.0	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	40.0
Purified water, ml	72.0	144.0	216.0	288.0	360.0	432.0	504.0	576.0	648.0	720.0	792.0	960.0

## Test procedure

**abia Progesterone** for the quantitative determination of progesterone concentration in human serum

---

- 1 Take the required number of coated strips. Place the unused strips back into the bag; reseal the foil-lined package in plastic zip-lock bag. Do not remove desiccant.

---

- 2 Analyze each calibrator, control serum, and sample duplicate.  
Add 25  $\mu$ l of calibrators 0 - 5 into appropriate wells.  
Add 25  $\mu$ l of control serum into appropriate wells.  
Add 25  $\mu$ l of samples to be tested in rest of the wells.  
The total time should not exceed 10 min.

---

- 3 Add 100  $\mu$ l conjugate into each well.  
Mix the contents of the wells for 30 seconds by careful tapping on the edge of the plate, then cover the plate with protective film.

---

- 4 Incubate for 90 minutes at room temperature 20–25 °C.

---

- 5 Remove the protective film slowly and carefully to prevent splashes. Aspirate the contents of all wells into a container for biohazardous waste (containing disinfectant).  
Add not less than 300  $\mu$ l of working washing solution into each well and aspirate. Perform this procedure 5 times. Use double aspiration in the final step where possible.

---

- 6 Add 100  $\mu$ l of TMB/substrate solution to all the wells. Keep the plates in a dark place for  $25 \pm 5$  minutes at 20–25 °C.

---

- 7 Add 150  $\mu$ l of stopping reagent into each well. Mix gently for 5–10 sec.

---

- 8 Read the optical density at 450 nm using a plate reader within 20 minutes after stopping reaction.

---

## Automated analyzers

Validated protocols for automated analyzers can be obtained from your representative. For the instrumentation without established validated protocol follow section “Test procedure” and ensure all requirements described in section “Precautions” are followed. All protocols for automated analyzers must be fully validated prior usage.

## Calculation and interpretation of the results

### Assay validation

Results of an assay are valid if the following criteria for the controls are met.

The absorbance (OD) of control serum should be within established range.  
The absorbance (OD) of calibrator 0 should be greater than 1.300.

### Calculation procedure

- 1 Calculate the mean optical density of each calibrator duplicate.
- 2 Calculate the mean optical density of each sample duplicate.
- 3 Draw a calibration curve on linear graph paper with the mean optical densities on the Y-axis and the calibrator concentrations on the X-axis.
- 4 Read the values of the unknowns directly off the calibration curve.  
If immunoassay software is being used, a 4-parameter curve is recommended.

Example	OD 1	OD 2	Mean OD	Value, nmol/l
Calibrator 0	2.496	2.464	2.480	0.00
Calibrator 1	2.193	2.201	2.197	1.50
Calibrator 2	1.396	1.367	1.382	10.00
Calibrator 3	0.801	0.789	0.795	25.00
Calibrator 4	0.508	0.497	0.503	50.00
Calibrator 5	0.297	0.295	0.296	100.00
Sample	0.780	0.785	0.783	26.20

This data is for illustration only and should **not be used** to calculate of samples. Each user should obtain his or her own data and standard curve.



## Performance characteristics

### Analytical sensitivity

The analytical sensitivity (limit of detection) was calculated by determining the variability of the calibrator 0 based on 10 replicate analyses less 2 x SD. Limit of detection defined at 0.5 nmol/l.

### Specificity

	Cross reactivity, %
Cortisone	0.041
Corticosterone	0.10
Testosterone	0.01
Estradiol (E2)	0.004
Dihydroprogesterone	0.041
Androstendione	0.10
Cortisol	0.005
17-OH-Progesterone	2.90
Pregнадione	10.00

### Precision

	Mean, nmol/l	SD	CV, %
Intra-assay, sample 1	25.98	1.130	4.37
Intra-assay, sample 2	4.75	0.300	6.26
Inter-assay, sample 1	26.50	1.600	6.10
Inter-assay, sample 2	5.10	0.300	6.80

### Accuracy

The assay was compared with a chemiluminescent microparticle immunoassay as a reference test. The total number of specimens was 380. The values ranged from 0.5 to 100 nmol/l. The least square regression equation and the correlation coefficient were computed for abia Progesterone in comparison with the reference method.

The least square regression analysis was  $y = 0.9689(x) + 4.02$  with correlation coefficient of 0.96.

### Expected normal value

	Range, nmol/l	
Males	0.50	5.20
Females, follicular phase	0.50	6.50
Females, luteal phase	8.00	87.00
Pregnant women, 1st trimester (1 - 12 weeks)	10.00	182.00
Pregnant women, 2nd trimester (13 - 24 weeks)	60.00	332.00
Pregnant women, 3rd trimester (25 - 40 weeks)	185.00	960.00

Normal value ranges may vary slightly among different laboratories. It is strongly recommended that each laboratory should determine its own range of expected normal values.













## Limitations of test

- the assay was validated only for the determination of progesterone in human serum
- the results obtained with this assay should never be used as the sole basis for clinical diagnosis. Any laboratory result is only a part of the total clinical picture of the patient
- the occurrence of heterophilic antibodies in patients regularly exposed to animals or animal products has to potential of causing interferences in immunological tests

## References

1. Buster J.E., R.J. Chang, D.L. Preston, et al: Interrelationships of circulating maternal steroids; progesterone, 16 -hydroxyprogesterone, 17 -hydroxyprogesterone, 20 -dihydroprogesterone, gamma-5-pregnenolone, gamma-5-pregnenolone-sulfate, gamma-5-pregnenolone-sulfate and 17-hydroxy gamma-5-pregnenolone, *J. Clin. Endocrinol. Metab.* 48:133, 1979.
2. Matthews C.P., et al.: *Obstet. Gynecol.*, 68:390, 1986.
3. Check, J.H., et al, Falsely elevated steroidal assay levels related to heterophile antibodies against various animal species. *Gynecol Obstet Invest* 40:139-140, 1995.

## Key to symbols used

	Manufacturer
	For in vitro diagnostic use
	Catalogue number
	Batch code
 YYYY-MM-DD	Expiry date
 2°C — 8°C	Storage temperature limitation
	Do not use if package is damaged
	Do not reuse
 $\Sigma_n$	Sufficient for [n] tests
	Consult Instructions for use
	Caution, consult documents
	Changes highlighted

## Hazard and precautionary statements for certain kit components

### Stopping reagent



Warning

H315	Causes skin irritation.
H319	Causes serious eye irritation.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### Conjugate, calibrators 0 - 5, control serum



Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P333 + P313	If skin irritation or rash occurs: Get medical advice/attention.

### Attention!

For complete precautionary statements and detailed information see safety data sheets (SDS).



**AB Diagnostic Systems GmbH**  
Sportfliegerstraße 4  
12487 Berlin  
Germany

☎ +49 30 208 987 160  
☎ +49 30 208 987 199  
✉ [info@ab-ds.de](mailto:info@ab-ds.de)  
[www.ab-ds.de](http://www.ab-ds.de)

